



## Extraction of RNA from Cells Cultured in Alvetex® Scaffold in 3D

### Introduction

The following protocol outlines the use of a commercial extraction kit (Qiagen RNeasy® kit, 74104) for the isolation of total RNA from cells cultured in Alvetex Scaffold. Example data was obtained using this protocol to extract RNA from HepG2 hepatocytes cultured in

Alvetex Scaffold for 7 days in 6-well inserts (AVP004-3) in Well Insert Holder in Deep Petri Dish (AVP015) format.

### Method

1. Wash Alvetex Scaffold culture by gentle immersion in PBS using flat-ended forceps and transfer to a clean 12-well plate.
2. Lyse cells by adding 600 µl Qiagen RNeasy® kit lysis buffer RLT per well.
3. Place on a rotating platform (100 rpm) for 10 minutes at room temperature.
4. Homogenise the lysate by passage up and down through a 20-gauge needle 10 times using a sterile plastic syringe.
5. Add an equal volume of 70 % ethanol to the homogenised lysate and pipette up and down 10 times to mix.
6. Transfer sample (up to 700 µl at a time) to an RNeasy® spin column in a 2 ml collection tube. Close cap and microcentrifuge for 15 seconds at 8000 x g. Discard flow-through. Repeat for remaining volume of lysate. Optional: On-column DNase digestion can be performed at this point (steps 6 a-d). Alternatively proceed to step 7.
  - a. Add 350 µl Buffer RW1 to the column. Close cap, and microcentrifuge for 15 seconds at 8000 x g to wash the membrane. Discard flow-through.
  - b. Prepare DNase 1 reagent by adding 10 µl DNase 1 stock solution (prepared according to manufacturer's instructions) to 70 µl Buffer RDD. Mix by gentle pipetting or gentle tube inversion.
  - c. Add the 80 µl diluted DNase 1 reagent to the column membrane and incubate at room temperature for 15 minutes.
  - d. Add 350 µl Buffer RW1 to the column. Close cap and microcentrifuge for 15 seconds at 8000 x g. Discard flow-through. Proceed to step 8.
7. Add 700 µl Buffer RW1 to the column. Close cap and microcentrifuge for 15 seconds at 8000 x g, reusing the collection tube. Discard flow-through.
8. Add 500 µl Buffer RPE (appropriately diluted with ethanol according to manufacturer's instructions) to the column. Close cap and microcentrifuge for 15 seconds at 8000 x g. Discard flow-through.
9. Add a further 500 µl Buffer RPE to the column. Close cap and microcentrifuge for 2 minutes at 8000 x g. Discard flow-through.
10. Transfer column to a clean 2 ml collection tube. Close cap and microcentrifuge at full speed



## Extraction of RNA from Cells Cultured in Alvetex® Scaffold in 3D

for 1 minute.

- Transfer column to a clean 1.5 ml collection tube and add 50 µl RNase-free water onto the column membrane. Close cap and leave at room temperature for 5 minutes. Microcentrifuge for 1 minute at 8000 x g to elute RNA.

### Example: Extraction of Total RNA from HepG2 Hepatocytes Cultured in Alvetex Scaffold in 3D

#### Cell Culture details

HepG2 cells (ATCC, HB-8065) were routinely maintained in T-75 flasks. HepG2 complete media consisted of: MEM media (Gibco, 21090) supplemented with 10 % v/v FBS, 2 mM L-glutamine and 100 U/ml Penicillin/Streptomycin. Cells were seeded onto Alvetex Scaffold discs in 6-well inserts (AVP004-3) in Well Insert Holders in Deep Petri dishes (AVP015), at a density of  $1 \times 10^6$  cells in 150 µl media suspension per insert. After settling for 3 hours in an incubator (5 % CO<sub>2</sub>, 37 °C) complete media was carefully added (70 ml per Petri dish). Cultures were maintained for 7 days, with a single media change on day 5. After 7 days, well inserts were dismantled and cultures were processed according to the protocol described above.


#### RNA analysis results

RNA samples were assessed for quantity and quality by spectrophotometric (Table 1) and agarose gel (Figure 1) analyses.

	Sample 1	Sample 2
RNA yield (µg/ml)	185.1 (±0.044)	182.7 (±0.470)
A <sub>260</sub> /A <sub>280</sub>	2.123 (±0.012)	2.117 (±0.006)

**TABLE 1.** RNA yield and purity analysis. RNA samples from duplicate HepG2 Alvetex Scaffold cultures were diluted 5-fold in 10 mM Tris-HCl, pH 7.0. Values shown represent the mean of triplicate readings per sample (±SD), taken on a NanoDrop 1000 spectrophotometer (Thermo Scientific).

## Extraction of RNA from Cells Cultured in Alvetex® Scaffold in 3D



**FIGURE 1.** Agarose gel RNA profile from duplicate HepG2 Alvetex Scaffold 3D cultures (Samples 1 and 2). 28S and 18S rRNA bands are clearly visible confirming RNA integrity (1% agarose gel; visualisation via ethidium bromide staining; 1 kb molecular weight ladder: Promega, G571A).

For research use only. Not for use in diagnostic procedures.



## Extraction of RNA from Cells Cultured in Alvetex® Scaffold in 3D

### Technical Support

To obtain the latest technical support information, please visit [www.reinnervate.com](http://www.reinnervate.com). You can have 24-hour access to the following information:

- Access to telephone and fax information for order support
- Search through the Reinnervate protocol video gallery
- Search product spec sheets, protocols, applications notes and product support docs
- Obtain information on customer training and workshops
- Submit a question directly to Technical Support at [tech.support@reinnervate.com](mailto:tech.support@reinnervate.com)

---

### Ordering Information

For the latest products and services available through Reinnervate, please visit [www.reinnervate.com](http://www.reinnervate.com).

#### Reinnervate

NETPark Incubator, Thomas Wright Way,  
Sedgefield, Co. Durham, TS21 3FD  
UNITED KINGDOM

T: +44 (0) 1740 625 266

F: +44 (0) 1740 625 251

E: [info@reinnervate.com](mailto:info@reinnervate.com)

---

### Trademarks

Unless otherwise noted, all trademarks are the property of Reinnervate Ltd.  
©2015 Reinnervate Ltd. All rights reserved.