

## Cost-Effective, Scalable, Lab-Scale iPS Spheroid Cultures



### Overview


Growth of iPS cells in spheroid suspension culture (3D) closely resemble embryoid body structures naturally formed by embryonic stem cells. This 3D growth format offers improved yield, viability and high efficiency for the expansion of human iPS stem cells and differentiation in a cost-effective format.

### The ABLE 30 mL Disposable Bioreactor

A delta-wing-shaped impeller with a magnet on each blade, located within the bioreactor, provides low-shear agitation by laminar flow, encouraging the formation and growth of uniform spheroid cell clusters.

The bioreactor allows for the flexibility of using seed stock from single-cell suspensions of iPS stem cells, passaged colony fragments, or mini-spheroids generated from low-adherence cell culture dishes or multi-well plates.

Bioreactor vessels of other sizes are available: 5 mL, 100 mL. The 5 mL size requires a different Stirrer Platform. A 500 mL size designed for a stand-alone bioreactor system, is also available.




**Human iPS cells (cell line 1231A3) grown in StemFit™ medium** demonstrating the consistency of spheroid sizes after 4 days cultivation in the ABLE 3D Disposable Bioreactor at 40 rpm.

### Benefits of Spheroid Culture

The ABLE 30 mL Disposable Bioreactor provides a combination of efficiency, embryoid body formation, cell differentiation, and cost savings. The Bioreactor comes ready-to-use for non-adherent cell growth, eliminating the need for high-priced, high performance extracellular matrix (ECM) proteins for coating plasticware. One bioreactor enables expansion for up to  $5 \times 10^7$  cells; approximately equivalent to cell yield from ten 10 cm culture dishes or ten 6-well plates. Ten dishes have approximately 550 cm<sup>2</sup> of combined surface area.

### Bioreactor Features

- Bioreactor capacity is 30 mL, enabling 5 to 10-fold cell expansion (up to  $5 \times 10^7$  cells)
- Typical growth is 200-300  $\mu$ m spheroids
- Bioreactor is made of high density polycarbonate for compatibility with iPS stem cell cultivation
- Vessel interior is surface treated for bio-compatibility
- Polypropylene screw cap 0.2  $\mu$ m nylon filter allows for passive gas exchange



### Single cell seeding provides for exponential growth

Human iPS cells (1231A3) were maintained on iMatrix-511 (Cat. No. NP891-011) – coated plates and grown in StemFit® medium (Cat. No. AKN02). Cells were harvested and dissociated into single cells using TrypLE Select (ThermoFisher), washed and counted. The single-cells were then seeded at 105 cells/mL in StemFit medium supplemented with 10  $\mu$ M Y27632 (Cat. No. 04-0012) and transferred to the ABLE 30 mL Disposable Bioreactor with constant spinner agitation at 55 rpm. Cells were harvested at day 2 and day 4, and spheroids dispersed by TrypLE Select, stained with trypan blue and counted.

ABLE company name and logo and are the property of ABLE Biotech Corp., Japan. REPROCELL name and logo are the property of REPROCELL Inc. © 2019 REPROCELL, Inc. All rights reserved.

# ABLE 3D Magnetic Stir and Disposable Bioreactor System for Stem Cell Culture



## Able Bioreactor Sizes



	5 mL	30 mL	100 mL	500 mL
Vessel Size:	5 mL	30 mL	100 mL	500 mL
Typical cell density:	10 <sup>6</sup> cells	10 <sup>7</sup> cells	10 <sup>8</sup> cells	10 <sup>9</sup> cells
Cat. No.:	ABBWVS05A	ABBWVS03A	ABBWVS10A	ABBWVS50A



Controller and Motor



Stirrer Platform

## ABLE System Configuration

The ABLE 3D System has only two components: a six-channel magnetic stirrer platform (bottom, 30 mL size shown), and an electrical controller and motor unit (top). The Stirrer Platform is placed inside a CO<sub>2</sub> cell culture incubator allowing the temperature, humidity and gas controls to be passively managed. The controller unit attaches magnetically to the outside of the incubator, and the motor is connected to the Stirrer Platform by a cable.

Cat. No.	Description	Size
ABBWDW-1013	ABLE Bioreactor Magnetic Stir System Controller and Motor	1 each
ABBWBP05N0S-6	ABLE Bioreactor Magnetic Stir System Base 5 mL	1 each
ABBWBP03N0S-6	ABLE Bioreactor Magnetic Stir System Base 30 mL and 100 mL	1 each
ABBWVS05A	ABLE 5mL Disposable Bioreactor	Pack of 6
ABBWVS03A-6	ABLE 30mL Disposable Bioreactor	Pack of 6
ABBWVS10A	ABLE 100 mL Disposable Bioreactor	Pack of 4
ABBWVS50A	ABLE 500 mL Disposable Bioreactor	1 unit

## References

- Matsuura K. Fabrication of mouse embryonic stem cell-derived layered cardiac cell sheets using a bioreactor culture system. *PLoS One* 7: e52176 (2012).
- Matsuura K. Creation of human cardiac cell sheets using pluripotent stem cells. *Biochem Biophys Res Commun* 425: 321-327 (2012).

### REPROCELL Europe Ltd

Thomson Pavillion, Todd Campus  
West of Scotland Science Park  
Acre Road  
Glasgow, G20 0XA  
UK

T: +44 (0)141 465 3460  
E: info-emea@reprocell.com

### REPROCELL USA Inc

9000 Virginia Manor Road  
Suite 207  
Beltsville, MD 20705  
USA

T: +1 301 470 3362  
E: info-us@reprocell.com

### REPROCELL India Ltd

3-1-135/1A, CNR Complex  
Mallapur  
Hyderabad 500 076  
Telangana  
India

T: +91-40-27178178  
E: Bhargavi.Gurram@reprocell.com

### REPROCELL Inc (Japan)

MetLife Shin-yokohama 381, Bldg. 9F  
3-8-11, Shin-yokohama  
Kohoku-ku, Yokohama  
Kanagawa 222-0033  
Japan

T: +81 45 475 3887  
E: info-asia@reprocell.com