

General Protocol for ICC Staining of Live Cells

Overview

The following procedure is a general protocol for immunocytochemistry (ICC) of live cells in culture. All steps should be performed in a sterile working environment using fresh culture medium for the continuation of the cell culture. Only open the antibody in the biological safety cabinet to prevent possible contamination.

Required Materials

- Cell culture medium

ICC Procedure

1. Culture the cells at 37°C and 5% CO₂ to the desired confluency.
2. Dilute the primary antibody to a final concentration of 2.5 to 5 µg/mL in fresh cell culture medium.
3. Aspirate the medium from the cells and add the diluted primary antibody to the live cells.
4. Incubate for 30 minutes at 37°C and 5% CO₂.
5. Aspirate the primary antibody and wash the cells gently 2 times with cell culture medium.
6. Add fresh cell culture medium and examine under a fluorescent microscope using the appropriate filters.

Note: Do not allow cultures to be kept outside of the incubator for extended periods of time to avoid drastic temperature changes and the possibility of contamination.

7. Return the plate to the incubator to continue culture at 37°C and 5% CO₂.

REPROCELL USA Inc

9000 Virginia Manor Road
Suite 207
Beltsville, MD 20705
USA

Tel: +1 301 470 3362

Email: info-us@reprocell.com

REPROCELL Europe Ltd

Thomson Pavilion, Todd Campus
West of Scotland Science Park
Acre Road Glasgow
G20 0XA
UK

Tel: +44 (0)141 465 3460

Email: info-emea@reprocell.com

REPROCELL India Ltd

3-1-135/1A, CNR Complex
Mallapur
Hyderabad 500 076
Telangana
India

Tel: +44 (0)141 465 3460

Email: Bhargavi.Gurram@reprocell.com

REPROCELL Inc (Japan)

MetLife Shin-yokohama 381, Bldg. 9F
3-8-11, Shin-yokohama Kohoku-ku,
Yokohama
Kanagawa 222-0033
Japan

Tel: +81 45 475 3887

Email: info-asia@reprocell.com

www.reprocell.com



PR-STEMGENT-ICC-LIVE-D001